

## Product Datasheet

# Human Ubiquitin Conjugating Enzyme E2A (UBE2A) ELISA Kit (orb1206774)

<b>Description</b>	Human Ubiquitin Conjugating Enzyme E2A (UBE2A) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.312-20ng/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Uniprot ID</b>	<a href="#">P49459</a>
<b>Sensitivity</b>	0.081ng/mL
<b>Entrez</b>	<a href="#">7319</a>
<b>Expiration Date</b>	Please enquire.