

Product Datasheet

Human Tau-Protein Kinase 1 (tPK1) ELISA Kit (orb1206782)

Description	Human Tau-Protein Kinase 1 (tPK1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9H3S4
Sensitivity	0.065ng/mL
Entrez	27010
Expiration Date	Please enquire.