

## Product Datasheet

### Human Relaxin (RLN) ELISA Kit (orb1206840)

<b>Description</b>	Human Relaxin (RLN) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	31.2-2000pg/mL
<b>Note</b>	For research use only
<b>Application notes</b>	competitive inhibition
<b>Sample Types</b>	serum, plasma or other biological fluids.
<b>Uniprot ID</b>	<b>P04808</b>
<b>Sensitivity</b>	9.66pg/mL
<b>Entrez</b>	<b>6013</b>
<b>Expiration Date</b>	Please enquire.