

Product Datasheet

Human Myeloid Progenitor Inhibitory Factor 1 (MPIF1) ELISA Kit (orb1206931)

Description	Human Myeloid Progenitor Inhibitory Factor 1 (MPIF1) ELISA Kit
Reactivity	Human
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	P55773
Sensitivity	5.4pg/mL
Entrez	6368
Expiration Date	Please enquire.