

Product Datasheet

Human Immune Receptor Expressed On Myeloid Cells 1 (IREM1) ELISA Kit (orb1207001)

Description	Human Immune Receptor Expressed On Myeloid Cells 1 (IREM1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q8TDQ1
Sensitivity	0.062ng/mL
Entrez	146722
Expiration Date	Please enquire.