

Product Datasheet

Human Intercellular Adhesion Molecule 4 (ICAM4) ELISA Kit (orb1207007)

Description	Human Intercellular Adhesion Molecule 4 (ICAM4) ELISA Kit
Reactivity	Human
Range	0.781-50ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q14773
Sensitivity	0.35ng/mL
Entrez	2286
Expiration Date	Please enquire.