

Product Datasheet

Human Heat Shock 70kDa Protein 1A (HSPA1A) ELISA Kit (orb1207010)

Description	Human Heat Shock 70kDa Protein 1A (HSPA1A) ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	PODMV8
Sensitivity	0.54ng/mL
Entrez	3303
Expiration Date	Please enquire.