

Product Datasheet

Human Gonadotropin Releasing Hormone Receptor (GnRHR) ELISA Kit (orb1207040)

Description	Human Gonadotropin Releasing Hormone Receptor (GnRHR) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P30968
Sensitivity	0.07ng/mL
Entrez	2798
Expiration Date	Please enquire.