

## Product Datasheet

# Human Growth Hormone Inducible Transmembrane Protein (GHITM) ELISA Kit (orb1207045)

<b>Description</b>	Human Growth Hormone Inducible Transmembrane Protein (GHITM) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.312-20ng/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Uniprot ID</b>	<b>Q9H3K2</b>
<b>Sensitivity</b>	0.13ng/mL
<b>Entrez</b>	<b>27069</b>
<b>Expiration Date</b>	Please enquire.