

Product Datasheet

Human Cholecystokinin 8, Octapeptide (CCK8) ELISA Kit (orb1207163)

Description	Human Cholecystokinin 8, Octapeptide (CCK8) ELISA Kit
Reactivity	Human
Range	15.6-1000pg/mL
Note	For research use only
Application notes	competitive inhibition
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sensitivity	4.46pg/mL
Expiration Date	Please enquire.