

## Product Datasheet

# Human Arylalkylamine-N-Acetyltransferase (AANAT) ELISA Kit (orb1207260)

<b>Description</b>	Human Arylalkylamine-N-Acetyltransferase (AANAT) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	1.56-100U/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	tissue homogenates, cell lysates or other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q16613</a>
<b>Sensitivity</b>	0.6U/mL
<b>Entrez</b>	<a href="#">15</a>
<b>Expiration Date</b>	Please enquire.