

Product Datasheet

Human Alpha-1-Antitrypsin (a1AT) ELISA Kit (orb1207261)

Description	Human Alpha-1-Antitrypsin (a1AT) ELISA Kit
Reactivity	Human
Range	31.2-2000ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P01009
Sensitivity	11.37ng/mL
Entrez	5265
Expiration Date	Please enquire.