

## Product Datasheet

# Rat Platelet Derived Growth Factor Subunit B (PDGFB) ELISA Kit (orb1207266)

<b>Description</b>	Rat Platelet Derived Growth Factor Subunit B (PDGFB) ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	31.2-2000pg/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates or other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q05028</a>
<b>Sensitivity</b>	12.5pg/mL
<b>Entrez</b>	<a href="#">24628</a>
<b>Expiration Date</b>	Please enquire.