

Product Datasheet

Rat Endothelin 2 (EDN2) ELISA Kit (orb1207284)

| | |
|--------------------------|--|
| Description | Rat Endothelin 2 (EDN2) ELISA Kit |
| Reactivity | Rat |
| Range | 7.81-500pg/mL |
| Note | For research use only |
| Application notes | sandwich |
| Sample Types | serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| Uniprot ID | P23943 |
| Sensitivity | 2.96pg/mL |
| Entrez | 24324 |
| Expiration Date | Please enquire. |