

## Product Datasheet

# Porcine Transforming Growth Factor Beta 1 (TGFb1) ELISA Kit (orb1207295)

<b>Description</b>	Porcine Transforming Growth Factor Beta 1 (TGFb1) ELISA Kit
<b>Reactivity</b>	Porcine
<b>Range</b>	31.2-2000pg/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, platelet-poor plasma, tissue homogenates, cell culture supernates or other biological fluids.
<b>Uniprot ID</b>	<b>P07200</b>
<b>Sensitivity</b>	10.5pg/mL
<b>Entrez</b>	<b>397078</b>
<b>Expiration Date</b>	Please enquire.