

Product Datasheet

Porcine Interleukin 5 (IL5) ELISA Kit (orb1207303)

Description	Porcine Interleukin 5 (IL5) ELISA Kit
Reactivity	Porcine
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9MYM5
Sensitivity	6.5pg/mL
Entrez	397409
Expiration Date	Please enquire.