

Product Datasheet

Bovine Interleukin 1 Receptor Antagonist (IL1RA) ELISA Kit (orb1207310)

Description	Bovine Interleukin 1 Receptor Antagonist (IL1RA) ELISA Kit
Reactivity	Bovine
Range	31.2-2000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	O77482
Sensitivity	13.5pg/mL
Entrez	281860
Expiration Date	Please enquire.