

Product Datasheet

Porcine Interleukin 1 Beta (IL1b) ELISA Kit (orb1207311)

Description	Porcine Interleukin 1 Beta (IL1b) ELISA Kit
Reactivity	Porcine
Range	62.5-4000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q29082
Sensitivity	26.2pg/mL
Entrez	396565
Expiration Date	Please enquire.