

Product Datasheet

Porcine Interleukin 1 Alpha (IL1a) ELISA Kit (orb1207313)

Description	Porcine Interleukin 1 Alpha (IL1a) ELISA Kit
Reactivity	Porcine
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P18430
Sensitivity	5.6pg/mL
Entrez	397094
Expiration Date	Please enquire.