

Product Datasheet

Bovine Interleukin 1 Alpha (IL1a) ELISA Kit (orb1207314)

Description	Bovine Interleukin 1 Alpha (IL1a) ELISA Kit
Reactivity	Bovine
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P08831
Sensitivity	5.8pg/mL
Entrez	281250
Expiration Date	Please enquire.