

Product Datasheet

Bovine Interleukin 18 (IL18) ELISA Kit (orb1207316)

Description	Bovine Interleukin 18 (IL18) ELISA Kit
Reactivity	Bovine
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9TU73
Sensitivity	6.1pg/mL
Entrez	281249
Expiration Date	Please enquire.