

Product Datasheet

Porcine Interleukin 17 (IL17) ELISA Kit (orb1207317)

Description	Porcine Interleukin 17 (IL17) ELISA Kit
Reactivity	Porcine
Range	31.2-2000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q60I29
Sensitivity	12.5pg/mL
Entrez	449530
Expiration Date	Please enquire.