

Product Datasheet

Porcine Interleukin 13 (IL13) ELISA Kit (orb1207323)

| | |
|--------------------------|--|
| Description | Porcine Interleukin 13 (IL13) ELISA Kit |
| Reactivity | Porcine |
| Range | 15.6-1000pg/mL |
| Note | For research use only |
| Application notes | sandwich |
| Sample Types | serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| Uniprot ID | Q95J68 |
| Sensitivity | 5.8pg/mL |
| Entrez | 396721 |
| Expiration Date | Please enquire. |