

Product Datasheet

Porcine Insulin Like Growth Factor 2 (IGF2) ELISA Kit (orb1207338)

Description	Porcine Insulin Like Growth Factor 2 (IGF2) ELISA Kit
Reactivity	Porcine
Range	1.56-100ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P23695
Sensitivity	0.64ng/mL
Entrez	396916
Expiration Date	Please enquire.