

Product Datasheet

Porcine Interferon Beta (IFN β) ELISA Kit (orb1207344)

| | |
|--------------------------|------------------------------------------------------------------------------------------------------|
| Description | Porcine Interferon Beta (IFN β) ELISA Kit |
| Reactivity | Porcine |
| Range | 15.6-1000pg/mL |
| Note | For research use only |
| Application notes | sandwich |
| Sample Types | serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| Uniprot ID | Q29412 |
| Sensitivity | 6.2pg/mL |
| Entrez | 445459 |
| Expiration Date | Please enquire. |