

Product Datasheet

Porcine Interferon Alpha (IFNa) ELISA Kit (orb1207346)

Description	Porcine Interferon Alpha (IFNa) ELISA Kit
Reactivity	Porcine
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P49879
Sensitivity	6.4pg/mL
Entrez	397686
Expiration Date	Please enquire.