

Product Datasheet

Bovine Interferon Alpha (IFNa) ELISA Kit (orb1207347)

Description	Bovine Interferon Alpha (IFNa) ELISA Kit
Reactivity	Bovine
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P05007
Sensitivity	5.9pg/mL
Entrez	282257
Expiration Date	Please enquire.