

Product Datasheet

Porcine Fibroblast Growth Factor 1, Acidic (FGF1) ELISA Kit (orb1207354)

Description	Porcine Fibroblast Growth Factor 1, Acidic (FGF1) ELISA Kit
Reactivity	Porcine
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	F1RM38
Sensitivity	6.6pg/mL
Entrez	397497
Expiration Date	Please enquire.