

Product Datasheet

Bovine Endocrine Gland Derived Vascular Endothelial Growth Factor (EG-VEGF) ELISA Kit (orb1207356)

Description	Bovine Endocrine Gland Derived Vascular Endothelial Growth Factor (EG-VEGF) ELISA Kit
Reactivity	Bovine
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	F1MN36
Sensitivity	6.3pg/mL
Entrez	617599
Expiration Date	Please enquire.