

Product Datasheet

Porcine Connective Tissue Growth Factor (CTGF) ELISA Kit (orb1207359)

Description	Porcine Connective Tissue Growth Factor (CTGF) ELISA Kit
Reactivity	Porcine
Range	0.625-40ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	F1S3P6
Sensitivity	0.242ng/mL
Entrez	100152141
Expiration Date	Please enquire.