

Product Datasheet

Porcine Bone Morphogenetic Protein 4 (BMP4) ELISA Kit (orb1207362)

Description	Porcine Bone Morphogenetic Protein 4 (BMP4) ELISA Kit
Reactivity	Porcine
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	A7LJT9
Sensitivity	0.067ng/mL
Entrez	100113425
Expiration Date	Please enquire.