

## Product Datasheet

# Bovine Angiotensin II (AngII) ELISA Kit (orb1207368)

|                          |   |
|--------------------------|---|
| <b>Description</b>       | Bovine Angiotensin II (AngII) ELISA Kit   |
| <b>Reactivity</b>        | Bovine                                    |
| <b>Range</b>             | 31.2-2000pg/mL                            |
| <b>Note</b>              | For research use only                     |
| <b>Application notes</b> | competitive inhibition                    |
| <b>Sample Types</b>      | serum, plasma or other biological fluids. |
| <b>Sensitivity</b>       | 8.46pg/mL                                 |
| <b>Expiration Date</b>   | Please enquire.                           |