

## Product Datasheet

### Porcine Adiponectin (ADP) ELISA Kit (orb1207369)

<b>Description</b>	Porcine Adiponectin (ADP) ELISA Kit
<b>Reactivity</b>	Porcine
<b>Range</b>	6.25-400ng/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma or other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q7YRF8</a>
<b>Sensitivity</b>	2.51ng/mL
<b>Entrez</b>	<a href="#">397660</a>
<b>Expiration Date</b>	Please enquire.