

## Product Datasheet

# Mouse Thrombopoietin (TPO) ELISA Kit (orb1207398)

|                          |   |
|--------------------------|---|
| <b>Description</b>       | Mouse Thrombopoietin (TPO) ELISA Kit      |
| <b>Reactivity</b>        | Mouse                                     |
| <b>Range</b>             | 78.1-5000pg/mL                            |
| <b>Note</b>              | For research use only                     |
| <b>Application notes</b> | sandwich                                  |
| <b>Sample Types</b>      | serum, plasma or other biological fluids. |
| <b>Uniprot ID</b>        | <a href="#">P35419</a>                    |
| <b>Sensitivity</b>       | 32pg/mL                                   |
| <b>Entrez</b>            | <a href="#">22018</a>                     |
| <b>Expiration Date</b>   | Please enquire.                           |