

Product Datasheet

Mouse Toll Like Receptor 6 (TLR6) ELISA Kit (orb1207405)

| | |
|--------------------------|---|
| Description | Mouse Toll Like Receptor 6 (TLR6) ELISA Kit |
| Reactivity | Mouse |
| Range | 0.312-20ng/mL |
| Note | For research use only |
| Application notes | sandwich |
| Sample Types | tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| Uniprot ID | Q9EPW9 |
| Sensitivity | 0.125ng/mL |
| Entrez | 21899 |
| Expiration Date | Please enquire. |