

Product Datasheet

Mouse Toll Like Receptor 4 (TLR4) ELISA Kit (orb1207406)

Description	Mouse Toll Like Receptor 4 (TLR4) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9QUK6
Sensitivity	0.065ng/mL
Entrez	21898
Expiration Date	Please enquire.