

Product Datasheet

Mouse Parathyroid Hormone Related Protein (PTHrP) ELISA Kit (orb1207444)

Description	Mouse Parathyroid Hormone Related Protein (PTHrP) ELISA Kit
Reactivity	Mouse
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P22858
Sensitivity	5.6pg/mL
Entrez	19227
Expiration Date	Please enquire.