

Product Datasheet

Mouse Prostaglandin Endoperoxide Synthase 1 (PTGS1) ELISA Kit (orb1207447)

Description	Mouse Prostaglandin Endoperoxide Synthase 1 (PTGS1) ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P22437
Sensitivity	0.126ng/mL
Entrez	19224
Expiration Date	Please enquire.