

Product Datasheet

Mouse Protein Z (PROZ) ELISA Kit (orb1207449)

| | |
|--------------------------|-------------------------------------|
| Description | Mouse Protein Z (PROZ) ELISA Kit |
| Reactivity | Mouse |
| Range | 15.6-1000pg/mL |
| Note | For research use only |
| Application notes | sandwich |
| Sample Types | plasma. |
| Uniprot ID | Q9CQW3 |
| Sensitivity | 6.5pg/mL |
| Entrez | 66901 |
| Expiration Date | Please enquire. |