

## Product Datasheet

# Mouse Pyruvate Kinase, Muscle (PKM2) ELISA Kit (orb1207458)

<b>Description</b>	Mouse Pyruvate Kinase, Muscle (PKM2) ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	0.312-20ng/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma or other biological fluids.
<b>Uniprot ID</b>	<a href="#">P52480</a>
<b>Sensitivity</b>	0.144ng/mL
<b>Entrez</b>	<a href="#">18746</a>
<b>Expiration Date</b>	Please enquire.