

## Product Datasheet

# Mouse Programmed Cell Death Protein 1 Ligand 1 (PDCD1LG1) ELISA Kit (orb1207473)

<b>Description</b>	Mouse Programmed Cell Death Protein 1 Ligand 1 (PDCD1LG1) ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	0.156-10ng/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma, tissue homogenates or other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q9EP73</a>
<b>Sensitivity</b>	0.058ng/mL
<b>Entrez</b>	<a href="#">60533</a>
<b>Expiration Date</b>	Please enquire.