

Product Datasheet

Mouse Olfactomedin 4 (OLFM4) ELISA Kit (orb1207485)

Description	Mouse Olfactomedin 4 (OLFM4) ELISA Kit
Reactivity	Mouse
Range	31.2-2000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	D3YYD0
Sensitivity	13.5pg/mL
Entrez	380924
Expiration Date	Please enquire.