

Product Datasheet

Mouse Matrix Metalloproteinase 15 (MMP15) ELISA Kit (orb1207517)

Description	Mouse Matrix Metalloproteinase 15 (MMP15) ELISA Kit
Reactivity	Mouse
Range	78.1-5000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	054732
Sensitivity	28.4pg/mL
Entrez	17388
Expiration Date	Please enquire.