

Product Datasheet

Mouse High Mobility Group Protein 1 (HMG1) ELISA Kit (orb1207570)

| | |
|--------------------------|--|
| Description | Mouse High Mobility Group Protein 1 (HMG1) ELISA Kit |
| Reactivity | Mouse |
| Range | 62.5-4000pg/mL |
| Note | For research use only |
| Application notes | sandwich |
| Sample Types | serum, plasma or other biological fluids. |
| Uniprot ID | P63158 |
| Sensitivity | 23.73pg/mL |
| Entrez | 15289 |
| Expiration Date | Please enquire. |