

Product Datasheet

Mouse Homing Associated Cell Adhesion Molecule (HCAM) ELISA Kit (orb1207575)

Description	Mouse Homing Associated Cell Adhesion Molecule (HCAM) ELISA Kit
Reactivity	Mouse
Range	78.1-5000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P15379
Sensitivity	29pg/mL
Entrez	12505
Expiration Date	Please enquire.