

Product Datasheet

Mouse Discoidin Domain Receptor Family, Member 1 (DDR1) ELISA Kit (orb1207631)

Description	Mouse Discoidin Domain Receptor Family, Member 1 (DDR1) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q03146
Sensitivity	0.047ng/mL
Entrez	12305
Expiration Date	Please enquire.