

## Product Datasheet

# Mouse Cross Linked C-Telopeptide Of Type II Collagen (CTXII) ELISA Kit (orb1207640)

<b>Description</b>	Mouse Cross Linked C-Telopeptide Of Type II Collagen (CTXII) ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	15.6-1000pg/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma or other biological fluids.
<b>Sensitivity</b>	4.26pg/mL
<b>Expiration Date</b>	Please enquire.