

Product Datasheet

Mouse Beta-Thromboglobulin (bTG) ELISA Kit (orb1207679)

Description	Mouse Beta-Thromboglobulin (bTG) ELISA Kit
Reactivity	Mouse
Range	78.1-5000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P62325
Sensitivity	34pg/mL
Entrez	12226
Expiration Date	Please enquire.