

Product Datasheet

Mouse Amyloid Beta Peptide 1-42 (Ab1-42) ELISA Kit (orb1207718)

| | |
|--------------------------|---|
| Description | Mouse Amyloid Beta Peptide 1-42 (Ab1-42) ELISA Kit |
| Reactivity | Mouse |
| Range | 15.6-1000pg/mL |
| Note | For research use only |
| Application notes | competitive inhibition |
| Sample Types | serum, plasma or other biological fluids. |
| Sensitivity | 5.48pg/mL |
| Expiration Date | Please enquire. |