

## Product Datasheet

### General 1-Desamino 8D Arginine Vasopressin (DDAVP) ELISA Kit (orb1207727)

<b>Description</b>	General 1-Desamino 8D Arginine Vasopressin (DDAVP) ELISA Kit
<b>Reactivity</b>	Other
<b>Range</b>	0.156-10ng/mL
<b>Note</b>	For research use only
<b>Application notes</b>	competitive inhibition
<b>Sample Types</b>	serum, plasma or other biological fluids.
<b>Sensitivity</b>	0.039ng/mL
<b>Expiration Date</b>	Please enquire.